# Maximum Performance in a Minimum of Space.

## MemoRail Modbus A1405 N

Compact multi-parameter measuring device for Memosens and digital sensors.



# Compact and Inexpensive.

MemoRail Modbus is a space-saving transmitter with Modbus interface in a 17.5 mm modular housing. Available as a 1- or 2-channel version for operation with contactless Memosens sensors as well as the SE 680 digital conductivity sensor and the SE 740 optical oxygen sensor.



## **Versatile Functionality**

The compact multi-parameter analyzer supports a variety of process variables:

- pH
- ORF
- Contacting conductivity
- Inductive conductivity
- Amperometric oxygen
- Optical oxygen
- Temperature

Red and green LEDs report the operating and sensor states. 1- and 2-channel versions are available to optimally suit your application.

### **Facts and Features**

- Slim modular housing with 17.5 mm width
- Memosens communication
- With pH/ORP combination sensors, measurement of up to four process values simultaneously plus temperature
- Modbus RTU protocol with standard RS-485 interface
- DTM for PACTware™ FDT or UnityPro FDT
- Up to 32 devices in parallel on a Modbus master
- Power supply 24 V DC
- 3-year warranty



## Easy to Use

For quick Modbus configuration, the network address is directly set on the MemoRail Modbus devices using the DIP switches on the front panel. 24 V DC is supplied via plug-in terminals or bus connectors.

DTM für PACTware™ FDT oder UnityPro FDT for easy handling of offline and online configurations. Available on www.knick.de

### **Memosens and Digital Sensors**

The use of digital sensors and Memosens sensors with contactless Memosens technology ensures maximum reliability and availability of the measuring point.

MemoRail Modbus is ready to use immediately after a Memosens or digital sensor is connected.

### **On-site Visualization**

Mobile tool for the contactless, visualization of Memosens measuring points without a display. Thanks to inductive data transmission, neither the process nor the measurement are interrupted. Easy and flexible evaluation of the on-site measuring point with the Portavo 907 and Portavo 908 multiparameter portables.

## **Universal Use**

MemoRail Modbus is an attractively priced solution in process analytics and optimally suited for use in numerous industries:

- Pharmaceutical industry, biotechnology
  - Upstream/downstream process
  - CIP/SIP systems
- Food & beverage
  - Process monitoring
  - CIP/SIP systems
- Water treatment plants/equipment
- Power generation, etc.

Because of their narrow modular housing, these DIN rail devices are ideal for installations where space is limited — in fermentation plants and in control cabinets, for example.





## Memosens

Memosens transfers both energy and data without contact between electrochemical sensors and analyzers. The integrated intelligence enables storing and analyzing process-related and sensor-relevant data directly in the sensors. These can be pre-calibrated under optimal laboratory conditions and quickly and easily installed at the measuring point or put in to replace sensors that have been "used up."

## **High Availability**

Through predictive maintenance, pre-calibrated digital Memosens sensors have a longer service life than analog sensors. The result is a significant reduction in maintenance costs and process downtime.

## **EMI Protection**

Perfect galvanic isolation ensures interference-free measurements, even under extreme potential conditions.

- Easy handling even under harsh conditions
- Up to 100 m cable length
- Resistant to moisture, dirt, corrosion, salt bridges, and interference potentials



# Wide Range of Sensors



MemoRail Modbus A1405 N can be combined with a large number of sensors, such as (from left to right): pH/ORP, conductivity, oxygen (optical), inductive conductivity, amperometric oxygen

## pH/ORP

Memosens sensors for measuring pH, ORP, and temperature. Different pH glasses and IsFET, reference systems, designs, and lengths make them perfectly adaptable to miscellaneous process requirements.

## Cond

2-electrode sensors with Memosens technology for measuring very low to medium conductivities.

Digital inductive conductivity sensors with extremely wide measuring range up to the highest conductivities. Range of applications from ultrapure water to concentration determinations.

## Оху

Amperometric Memosens sensors for measuring very low oxygen values up to pure oxygen, dissolved in water or gaseous.

Flow-independent optical oxygen sensor with fast response time. For measurements in hygienic areas; steam-sterilizable, autoclavable, and CIP-resistant.



# Convenient Calibration





## **Easy Sensor Calibration via Modbus**

MemoRail Modbus A1405 N offers the option to calibrate the sensors directly on site using the Modbus. A variety of calibration procedures is available for the various parameters.

Easy sensor calibration via the process control system with the free of charge DTM for PACTware™ FDT or UnityPro FDT. Available on www.knick.de.

## **Mobile Sensor Calibration**

A mobile alternative for offline calibration and adjustment of Memosens sensors is provided by the Portavo 908 Multi analyzer from Knick, which can also be used directly on site. The integrated temperature detector of the Memosens sensor can be very easily calibrated with the Portavo 908 Multi.

The Portavo 908 Multi portable multi-parameter analyzer for measuring pH, ORP, conductivity, or oxygen is available as a GLP-compliant version with printer interface for use in the pharmaceutical and biotechnology industries.







## **Sensor Calibration via MemoSuite**

The flexible and intuitive MemoSuite software tool enables easy calibration of Memosens sensors in the lab. On-site calibration under adverse ambient conditions is no longer necessary. The only thing required there is quick and uncomplicated replacement of the used sensors with pre-calibrated sensors.

Extensive buffer management is offered by the convenient compilation of individual buffer sets from an extensive library. It is also a simple matter to enter special, user-specific buffer tables incl. temperature values.

To meet a variety of application-specific requirements, MemoSuite is available in different versions:

- **MemoSuite Basic** for calibration of Memosens sensors.
- MemoSuite Advanced for calibration, diagnostics and database documentation of the sensors.
  Up to 10 sensors can be calibrated simultaneously. The database complies with the requirements of GMP and FDA CFR 21 Part 11; the complete documentation can be output either as a calibration report or as a dataset in Excel format.

## Knick >



Interface Technology **Indicators Industrial Transmitters Portables Laboratory Meters** Sensors **Fittings** 

## Knick The Art of Measuring

Knick has been among the leading manufacturers of electronic measurement devices for process analysis for more than 70 years. Today, the company alignment is still focused on a high technical level and an intensely innovative orientation.

The current product range of Knick includes unique device series also for extremely difficult applications. As a system provider, Knick can guarantee expert analysis of your measurement requirements and provide qualified advice for complete measurement loops. www.knick.de

## Knick Elektronische Messgeräte **GmbH & Co. KG**

Beuckestrasse 22, 14163 Berlin, Germany Phone: +49 30 80191-0

Fax: +49 30 80191-200 info@knick.de · www.knick.de